

Liquid Single-Use Bags

High-performance fluid storage and transfer bags with optimized reliability, throughput, and quality assurance for critical bioprocessing for a variety of upstream and downstream processes

Liquid single-use bags are low-profile containers engineered to bioprocess fluids. Featuring high-integrity and flow-optimized TruFLO ports, liquid single-use bags offer reliable and high-throughput performance. The industry-proven, medical-grade Renolit 9101 multilayer polyethylene film meets the requirements of ISO and USP biocompatibility tests. The film also meets low-permeability and low-temperature requirements while offering high clarity.

Liquid single-use bags are available in standard sizes from 50 mL to 20 L in 2 and 3-port options, and can be customized up to 50 L with up to 4-ports in larger chamber sizes (2 L to 50 L). Custom face-ported 2D bags are also available in 20 L, 50 L, 100 L and 200 L.

Liquid single-use bags feature a variety of BPOG-compliant standard and custom configurations

Feature	Configurations
Volume	Standard: 50 mL to 20 L Custom: up to 200 L
Port	Standard Edge: 2- or 3-port (500 mL to 50 L) Custom Edge: 4-port (2 L to 50 L) Custom Face: 3-port (20 L, 50 L, 100 L, 200 L)
Port Size	Edge: 1/8 in, 1/4 in, or 3/8 in Face: 1/4 in, 3/8 in or 1/2 in
Tubing	Standard: Thermoplastic Elastomer (TPE) Custom: Platinum-Cured Silicone (PCS)



TruFLO ports are specially engineered to offer high-integrity and improved flow for critical bioprocessing applications

Benefits

- Unique TruFLO ports design facilitates optimized edge-seal integrity and improved flow rates
- Low-profile design ensures minimal product holdup to maximize product recovery
- Configurable with a wide variety of BPOG-compliant components

Typical Applications

- Buffer and cell culture media
- Bulk product collections and storage
- Chromatography and filtration buffer
- Fraction collection
- Product sampling and transport

Standard Configurations

Liquid single-use bags with TruFLO ports



2-Port Configurations

Size	Port Sizes	Dimensions (W x L)	Internal Surface Area	Connections	Smart Part #
50 mL	1/4 in	5.37 x 4.25 in	20 in ²	Luer (Body & Insert)	EP050M2PTL-G
100 mL	1/4 in	5.37 x 4.75 in	25 in ²	Luer (Body & Insert)	EP100M2PTL-G
250 mL	1/4 in	5.37 x 6.50 in	41 in ²	Luer (Body & Insert)	EP250M2PTL-G
500 mL	1/4 in	7.07 x 7.75 in	65 in ²	Luer (Body & Insert)	EP500M2PTL-G
1 L	1/4 in	7.07 x 10.31 in	104 in ²	Luer (Body & Insert)	EP001L2PTL-G
500 mL	3/8 in	7.07 x 7.75 in	65 in ²	MPCs (F&M)	EP500M2PTM-G
1 L	3/8 in	7.07 x 10.31 in	104 in ²	MPCs (F&M)	EP001L2PTM-G
2 L	3/8 in	7.07 x 15.75 in	173 in ²	MPCs (F&M)	EP002L2PTM-G
5 L	3/8 in	14.25 x 14.91 in	280 in ²	MPCs (F&M)	EP005L2PTM-G
10 L	3/8 in	16.75 x 18.88 in	440 in ²	MPCs (F&M)	EP010L2PTM-G
20 L	3/8 in	17.50 x 28.38 in	772 in ²	MPCs (F&M)	EP020L2PTM-G



3-Port Configurations

Size	Port Sizes	Dimensions (W x L)	Internal Surface Area	Connections*	Smart Part #
500 mL	3/8 in (2) & 1/4 in	7.0x x 7.75 in	65 in ²	MPCs (F&M), Injection Cap	EP500M3PTM-G
1 L	3/8 in (2) & 1/4 in	7.07 x 10.31 in	104 in ²	MPCs (F&M), Injection Cap	EP001L3PTM-G
2 L	3/8 in (2) & 1/4 in	7.07 x 15.75 in	173 in ²	MPCs (F&M), Injection Cap	EP002L3PTM-G
5 L	3/8 in (2) & 1/4 in	14.25 x 14.91 in	280 in ²	MPCs (F&M), Injection Cap	EP005L3PTM-G
10 L	3/8 in (2) & 1/4 in	16.75 x 18.88 in	440 in ²	MPCs (F&M), Injection Cap	EP010L3PTM-G
20 L	3/8 in (2) & 1/4 in	17.50 x 28.38 in	772 in ²	MPCs (F&M), Injection Cap	EP020L3PTM-G

* via 12-inch TPE tubing

Medical-Grade Renolit 9101 PE Film




All liquid single-use bags are manufactured with industry-proven, medical-grade Renolit 9101 film, a multilayer film composed of a high-purity biocompatible polyethylene (PE) contact layer with internal ethylene vinyl alcohol (EVOH) oxygen-barrier.

Physical Properties	
Property	Typical Value*
Film Thickness	0.325 mm
Clarity	97% (ASTM D-1003)
Tensile Strength at Break	13MPa (ASTM D-882)
Elongation at Break	350% (ASTM D-882)
Break at Cold Temperature	< -45°C (ISO 8570)
Water Vapor Transmission [†]	0.32 g/m ² /day (ASTM F-1249)
O ₂ Permeability [†]	<0.05 cm ³ /m ² /day/bar (ASTM D-3985)
CO ₂ Permeability [†]	<0.2 cm ³ /m ² /day/bar (ASTM F-2476)

* Transmission values for film gamma-irradiated with 50 KGy. Other are for film gamma-irradiated with 25 KGy.
[†] @ 23 °C, 100% RH.
[‡] @ 23 °C, 0% RH.

Renolit 9101 PE Film Layers




Quality, Regulatory, and Biocompatibility Properties	
Category	Property/Test*
Composition	<ul style="list-style-type: none"> High-purity polyethylene (PE) and ethylene vinyl alcohol (EVOH) Animal-Derived Ingredient (ADI) Free
Biocompatibility	<ul style="list-style-type: none"> ISO 10993-4, Hemolysis ISO 10993-5, Cytotoxicity ISO 10993-6, Implantation ISO 10993-10, Irritation and Sensitization ISO 10993-11, Acute System Toxicity USP <85>, Bacterial Endotoxins – LAL test USP <87>, Biological Reactivity in vitro USP <88>, Biological Reactivity in vivo, Class VI
Extractables/Leachables	<ul style="list-style-type: none"> USP <661.1>, Polyethylene Physicochemical Tests, Extractable Metals, Plastic Additives Ph. Eur. 3.1.5, Polyethylene with additives for containers for parenteral preparations and for ophthalmic preparations

Recommended Sterilization Method: Gamma

* Pharmacopoeia and Biocompatibility compliance test reports available upon request

Questions or Inquiries?

More information is available at www.ilcdoover.com, or by contacting us at customer_service@ilcdoover.com or simply by reaching out to your dedicated ILC Dover sales representative.



BEYOND BOUNDARIES®

ILC Dover is a world-leader in the innovative design and production of engineered flexible protective solutions for pharmaceutical and biopharmaceutical, flood protection, personal protection, bulk packaging, and aerospace industries. Our customers will attest to our relentless dedication to high value products, advanced technology, and responsive service, as our visionary flexible protective solutions have improved efficiency while safeguarding people, product, and infrastructure in hazardous conditions since 1947.

 **HEADQUARTERS**
One Moonwalker Road
Frederica, DE (USA) 19946
customer_service@ilcdover.com
+1.302.335.3911 x506 - Tel
+1 302.335.1320 - Fax

 **CUSTOMER SERVICE**
+1.800.631.9567 (US Toll Free)

 **WEBSITE**
www.ilcdover.com

