



Gasket & Sealing Solutions

Sentinel Process Systems, Inc. has partnered with leading manufacturers of high quality sealing products to provide a wide variety of sanitary and hygienic seals. All of the products we represent are proudly made in the USA under the highest quality standards. Sentinel has been a critical partner to the Biotech, Pharmaceutical and Food & Beverage industry since 1950 and we support their operations with the largest and most diverse inventory in the country.

- Broad range of Gaskets and O-Rings
- Fully engineered, tested and lot traceable materials
- Standard and custom sealing designs ensure the highest purity and performance
- Color coded Gaskets and O-Rings
- Made in the U.S.A. for quick lead times
- FDA and USP Class VI

Visit us at [sentinelprocess.com](https://www.sentinelprocess.com) or call 800.345.FLOW today!

Contact Sentinel Process Systems, Inc. to help with all your Process Component, Instrumentation and Single-Use needs.

Gasket & Sealing Solutions



Tri-Clamp Gaskets
EPDM, Silicone, PTFE



Envelope Gaskets
Teflon™ with EPDM or Viton Filler



Bevel Seat Gaskets
Buna, FKM/Viton, EPDM,
PTFE, Platinum Silicone



Metal Detectable Gaskets
PTFE Blended Material



Screen Gaskets
FKM, Buna, EPDM, PTFE, Platinum Silicone
with 304 or 316L Stainless Steel Screens



Sock Screen Gaskets
Buna, FKM/Viton, EPDM,
PTFE, Platinum Silicone



Orifice Plates



O-Rings
EPDM, Platinum Cured Silicone,
FKM/Viton, Buna, PTFE



Delivering Innovative Process Technologies



800.345.FLOW
sentinelprocess.com