



Class VI, platinum-cured silicone elastomer comprises AdvantaPure's peristaltic pump grade silicone tubing. Ideally suited for pharmaceutical, biotech, and laboratory environments, AdvantaSil APSPG is specially formulated to withstand repeated compression and release for consistent, dependable performance – far outlasts standard silicone tubing. AdvantaSil APSPG pump tubing has undergone extensive physical, chemical, and biological testing, and meets USP Class VI, FDA CFR 177.2600, ISO 10993 and European Pharmacopoeia 3.1.9 standards.

AdvantaSil APSPG is not intended for implantation, continuous steam applications or elevated pressure levels.



KEY FEATURES

- Special formulation promotes extended life in pumps
- Extremely flexible and resilient; hardness value of 50 Shore A
- Outstanding compression characteristics – snaps back to original form
- Smooth surfaces offer uniform fluid flow and resistance to bacteria build-up
- AdvantaSil APSPG tubing is USP Class VI approved
- Material meets FDA 21 CFR 177.2600
- Certified free of animal-derived ingredients
- Low volatile grade suitable for pharmaceutical and biopharm applications
- Maximum operating temperature: 400°F (204.4°C); brittle temperature: less than -100°F (-73.3°C)
- Sterilizable/autoclavable
- Translucent for visual contact with the flow
- Stocked in sizes from .031" through .500" I.D.
- Documented lot traceable with identification on bags
- Documented quality control
- Custom lengths, sizes, and packaging available
- AdvantaFlex® biopharmaceutical grade tubing also available – heat sealable, weldable and moldable



SILICONE
www.advantapure.com




3265 Sunset Lane, Hatboro, PA 19040
www.sentinelprocess.com
P: 800.345.FLOW F: 888.329.9669

Complete validation package available on request



AdvantaSil silicone tubing and hose products are quickly becoming recognized as superior choices for Single-Use processes worldwide. Engineers in process development and manufacturing have relied on AdvantaSil unreinforced tubing, reinforced hose, pump grade tubing and wire reinforced suction hose for years.



SPECIFICATIONS

Product Number	I.D. (in.)	I.D. (mm.)	Wall (in.)	Wall (mm.)	O.D. (in.)	O.D. (mm.)	Standard Available Length (ft.)
APSPG-0031-0156	.031	.79	.063	1.59	.156	3.97	25
APSPG-0063-0188	.063	1.59	.063	1.59	.188	4.76	25
APSPG-0125-0250	.125	3.18	.063	1.59	.250	6.35	50
APSPG-0188-0313	.188	4.76	.063	1.59	.313	7.94	50
APSPG-0250-0375	.250	6.35	.063	1.59	.375	9.52	50
APSPG-0250-0438	.250	6.35	.094	2.38	.438	11.11	50
APSPG-0250-0500	.250	6.35	.125	3.18	.500	12.70	50
APSPG-0313-0500	.313	7.94	.094	2.38	.500	12.70	50
APSPG-0375-0563	.375	9.52	.094	2.38	.563	14.29	50
APSPG-0375-0625	.375	9.52	.125	3.18	.625	15.88	50
APSPG-0500-0750	.500	12.70	.125	3.18	.750	19.05	50

Sold by standard coil length only. Add length suffix code to product number when ordering—see the Coil Length Legend at right for length suffix codes. Example: 25 ft. of .031 in. I.D. x .156 in. O.D. tubing is product number APSPG-0031-0156P.
Coils are supplied double bagged in heat-sealed polybags and bulk packed. Contact your AdvantaPure Sales Representative for other packaging options.

Coil Length Legend:

P = 25 feet
L = 50 feet

SILICONE

Purity in Fluid Flow Systems®

AdvantaSil APSPG is part of a complete line of products for the Pharmaceutical, Biopharm, Chemical, and Food industries. Items available include . . .

- Unreinforced, 50 Shore A, platinum-cured silicone tubing
- Braid-reinforced, platinum-cured silicone hose
- Wire and mesh reinforced, mandrel-wrapped, platinum-cured silicone suction hose
- Stainless steel overbraided PTFE hose; rubber-covered FEP and EPDM hose
- Weldable and sealable AdvantaFlex® biopharmaceutical grade TPE tubing
- Sanitary fittings and complete hose assemblies including molded silicone ends
- Silicone injection molding capabilities – molded manifold assemblies & BioClosure® container closures



Address: NewAge Industries AdvantaPure®
145 James Way
Southampton, PA 18966 U.S.A.

Toll Free Phone: **888-755-4370**
Phone: **215-526-2151**
Toll Free Fax: **888-258-4293**
Fax: **215-526-2167**
Web Site: **www.advantapure.com**
E-mail: **sales@advantapure.com**



DISTRIBUTED BY:



3265 Sunset Lane, Hatboro, PA 19040
www.sentinelprocess.com
P: 800.345.FLOW F: 888.329.9669

©2017 NewAge® Industries, Inc. Specifications are subject to change without notice. AdvantaPure products are produced by NewAge® Industries, Inc. AdvantaSil™, AdvantaFlex®, BioClosure®, NewAge Industries AdvantaPure®, Purity in Fluid Flow Systems®, NewAge®, and the "N" logo are trademarks of NewAge® Industries, Inc.



WARRANTY: NewAge Industries AdvantaPure products are warranted to be free from defects in material and workmanship. Any product found to carry such defects will be replaced without charge. This warranty is limited only to the replacement of NewAge Industries AdvantaPure products, item for item, which after examination by NewAge Industries AdvantaPure are deemed defective. NewAge Industries AdvantaPure specifically disclaims any other liability. The data supplied is provided as a helpful guideline and is believed to be reliable; however, nothing stated shall constitute a guarantee, recommendation, or warranty for any application.

NOTICE: NewAge Industries, Inc. has been a privately owned corporation since 1954. Neither the corporation nor its principals are affiliated, in any way, with any other group or entity either representing themselves as or being referred to by the term NewAge.