

Notification of Pall Switching to Silicone Caps

Date: July 15, 2022

PCN#: CL-22-00475

To Our Valued Customers:

Entegris is committed to improving the quality, purity, and performance of our products. This notification is to inform you that Pall is switching to platinum cured silicone caps used on their filter offerings.

Change Description

Pall is switching to platinum cured silicone caps, and replacing PVC protective caps that are used on SUS component filters and Pall Capsule filters. There is no impact on the single-use system's form, fit, function, or product claims.

Reason for Change:

Pall was made aware that in very few instances, the tubing of some Allegro SUS showed some areas of discoloration. This has been attributed to the presence of the PVC caps. The replacement platinum cured silicone caps have qualified to show a PVC equivalent protection level.

Estimated Timing of Change:

This change is anticipated to be implemented the beginning of September 2022.

Entegris strives to provide you with the highest quality products and attempts to make changes with as little disruption to your process as possible. Therefore, a copy of this notification should be forwarded immediately to the personnel within your organization affected by this change.

Please communicate through your Entegris account manager, customer service or the change control administration group at pcn.change.control@entegris.com.

Thank you for your loyal business with Entegris.

Sincerely,
Chance Nguyen
Associate Product Manager – Life Sciences

Products Affected by PCN # CL-22-00475

Material	Material Description
FCH04516.01	FLTR MNFLD W/WASTE BG
FCH04834.01	Manifold, Sampling Assy BRX 50 ML
FCC03890.01	20 Liter, 3 Port with Filter
FCH03865.01	Filter Assembly with Pall Filter