



Single-Use Media Storage

Sentinel Process Systems, Inc. offers a wide range of Media Storage solutions used in single-use applications. Our fabrications are designed to meet complex demands for material integrity throughout the process. Available in standard designs or customized solutions to meet specific customer requirements with short lead times.

- Variety of storage designs to address every application
- Various materials including LDPE, EVOH, PTEG, PP, HDPE and PC
- Silicone and TPE tubing with numerous end connection options
- Gamma irradiated between 25 and 40kGy
- Stocking programs at no cost
- Assembled in ISO 7 cleanroom
- All product contact materials USP Class VI Tested
- Full traceability for all materials
- Certificates of analysis provided for each lot
- Extractables and leachables (E&L) data available for all contact materials

Visit us at [sentinelprocess.com](https://www.sentinelprocess.com) or call 800.345.FLOW today!

Contact Sentinel Process Systems, Inc. to help with all your Process Component, Instrumentation and Single-Use needs.


Single-Use Media Storage



2D & 3D Bioprocess Bags
50ml - 2,500L
1 - 5 Ports, Single, 2 and 3-Ply



Motion Bioreactor Bags
5L to 50L, 4, 5 and 7 Ports,
Single, 2 and 3-Ply



Cryopreservation Bags
2 or 3 Port, -321°F to 250°F



Bag Manifolds
Custom Sizes, Lengths
and Configurations



Tank Liners
50L to 1,500L Open or
Closed Top, 3-Ply



Totes



Bottles & Carboys
From 60ml to 40L



**Bottle Cap
Assemblies**



**Also Check Out Our
Tubing Solutions Flyer**



Delivering Innovative Process Technologies



800.345.FLOW
sentinelprocess.com